

RAYSHAPE

ShapeCure D mini

Complete and Even Curing

Professional Finish for Dental 3D Prints



360° Auto Rotation



Intuitive Operation



Low Maintenance



Validated Solution



ShapeCure D mini Curing Station delivers strong, stable 3D-printed dental applications — including crowns, bridges, dentures, surgical guides, splints, and more. Simple to run and reliably consistent, it ensures complete curing with multi-angle UV light. Get clinically accurate, ready-to-use parts right at your desk — effortlessly.

Specifications

Curing Volume	D160 mm × H135 mm	Wavelength	405 nm
Curing Power	90 W	Rated Power	120 W
Curing Setting	Customized duration	Control Chip	STM 32 ARM
Device Size	221 × 220 × 246 mm	Device Weight	3.1 kg
Package Size	285 × 285 × 320 mm	Package Weight	4.5 kg
Power Requirement	Adapter input: 100 ~ 240VAC, 50/60Hz, Max 1A; Rated input: DC 24V, 5A		

Contact Us

www.rayshape3d.com
sales.os@rayshape3d.com

Headquarters

Suzhou Laisai Intelligence Technology Co., Ltd.
No.3 Xinghua Parkway, Changshu, Jiangsu, China

EU Subsidiary

RAYSHAPE B.V.
William Boothlaan 258 Amstelveen The Netherlands

US Subsidiary

RAYSHAPE LLC
6605 Longshore Street, Suite 240, Dublin OH 43017, USA